Self-Cured Dental Universal Adhesive

Tokuyama Universal Bond II

 REF AD04

Instructions for Use

Manufacturer

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com/>

EC Rep.

Tokuyama Dental Italy S.r.l.

Via Chizzalunga, 1

36066 Sandrigo, Vicenza, Italy

Tel: +39 0444 659650

Fax: +39 0444 750345

Importer

Tokuyama Dental Deutschland GmbH

Am Landwehrbach 5

48341 Altenberge Germany

Tel: +49 2505 938513

Fax: +49 2505 938515

ENGLISH	3	ČESKY	12	DANSK	22	DEUTSCH	32
EESTI KEEL	43	FRANÇAIS	53	ITALIANO	63	LATVISKI	74
LIETUVIŲ K.	84	MAGYAR	94	NEDERLANDS	104	NORSK	114
POLSKI	124	SLOVENČINA	135	SUOMI	145	SVENSKA	155

MD	(EN) Medical Device	(CZ) Zdravotnický prostředek	(DK) Medicinsk udstyr
	(DE) Medizinprodukt	(EE) Meditsiiniseade	(FR) Dispositif médical
	(IT) Dispositivo medico	(LV) Mediciniska ierīce	(LT) Medicinos priemonė
	(HU) Orvostechnikai eszköz	(NL) Medisch hulpmiddel	(NO) Medisinsk utstyr
	(PL) Wyrób medyczny	(SK) Zdravotnícka pomôcka	(FI) Lääkinnällinen laite
	(SE) Medicinteknisk produkt		

This product is certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

1. Tokuyama Universal Bond II is a two-component self-cured dental adhesive system for both direct and indirect restorations that can be used with Self-etch, Selective-enamel-etch and Total-etch techniques. As a universal adhesive, Tokuyama Universal Bond II has been designed to be fully compatible with light-cured, self-cured and dual-cured composite materials.
2. Tokuyama Universal Bond II improves the bond strength of polymerizable resin material (adhesive resin cement, acrylic resin and composite resin) to indirect restorative materials such as glass-ceramics (porcelain), oxide-ceramics (zirconia and alumina), metals (precious and non-precious) and resin materials including inorganic filler.
3. Tokuyama Universal Bond II contains Phosphoric acid monomer, Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), 2-Hydroxyethyl methacrylate (HEMA), MTU-6 (thiouracil monomer), Silane coupling agent, Peroxide, Borate catalyst, Acetone, Ethanol and Purified water.

■ INDICATIONS

- Direct anterior and posterior restorations with light-curing, dual-curing, and self-curing composite materials
- Intraoral repair of composite restorations, porcelain fused to metal, metal, and all-ceramic restorations without an additional primer
- Cementation of indirect restorations and veneers when combined with light-cure, dual-cure, and self-curing resin cements
- Bonding of core build-ups made of core build-up materials
- Bonding of denture resin to metal base, clasp or attachment
- Repair of denture with metal base, clasp or attachment
- Bonding of opaque resin to a metal base in the fabrication of resin-faced crowns

■ CONTRAINDICATIONS

Tokuyama Universal Bond II contains methacrylic monomers, peroxide, organic solvents and acids. DO NOT use Tokuyama Universal Bond II for patients allergic to or hypersensitive to methacrylic monomers, related monomers, peroxide, organic solvents and acids.

■ PRECAUTIONS

- 1) DO NOT use Tokuyama Universal Bond II for any purpose other than those listed in these instructions. Use Tokuyama Universal Bond II only as directed herein.
- 2) Tokuyama Universal Bond II is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- 3) DO NOT use Tokuyama Universal Bond II if the safety seals are broken or appear to have been tampered with.
- 4) If Tokuyama Universal Bond II causes an allergic reaction or oversensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling Tokuyama Universal Bond II to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If Tokuyama Universal Bond II comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of Tokuyama Universal Bond II with eyes, mucosal membrane, skin and clothing.
 - If Tokuyama Universal Bond II comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If Tokuyama Universal Bond II comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the completing the restoration. Affected areas may whiten from protein coagulation; however, whitening should disappear within 24 hours. If whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
 - If Tokuyama Universal Bond II comes into contact with the skin or clothing, immediately saturate the area with an alcohol-soaked cotton swab or gauze and flush with water.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) Tokuyama Universal Bond II should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 8) To avoid the unintentional ingestion of Tokuyama Universal Bond II, do not leave it unsupervised within the reach of patients and children.
- 9) DO NOT expose Tokuyama Universal Bond II or its vapor to open flame because Tokuyama Universal Bond II is flammable.
- 10) To avoid cross infection or lowering bond strength, DO NOT reuse the disposable applicator

- and the disposable mixing well. Clean the mixing well thoroughly with alcohol after each use.
- 11) DO NOT use Tokuyama Universal Bond II directly to the cavity being in close proximity to the pulp. Pulp protection with glass ionomer lining or calcium hydroxide is recommended. DO NOT use EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of Tokuyama Universal Bond II.
 - 12) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
 - 13) DO NOT use mixing wells supplied by other manufacturers.
 - 14) Use the mixing well or disposable mixing well provided when dispensing and mixing of Tokuyama Universal Bond II. If Tokuyama Universal Bond II is not dispensed into the rubber, it may be difficult to remove the cured Tokuyama Universal Bond II from mixing well. The cured Tokuyama Universal Bond II can be removed by inverting the rubber.
 - 15) DO NOT mix Tokuyama Universal Bond II with other brands of primers or adhesives.
 - 16) To avoid splashing the liquid, the bottle of Tokuyama Universal Bond II should not be opened immediately after taking out of the refrigerator. Tokuyama Universal Bond II should be brought to room temperature before using. After removing from the refrigerator, allow the product to sit for 20 minutes or until it reaches room temperature.
 - 17) Do not use Tokuyama Universal Bond II with resin materials for which special adhesives or primers are specified.

■ PRECAUTIONS FOR MEDICAMENTS AND MATERIALS

- 1) Some materials and medicaments (hemostatic material) inhibit setting/adhesion of Tokuyama Universal Bond II.
DO NOT USE products that contain:
 - eugenol,
 - aluminum chloride,
 - ferric sulfate,
 - calcium sulfate,
 - aluminum sulfate,
 - diammine silver fluoride [molecular formula: $\text{Ag}(\text{NH}_3)_2\text{F}$]
- 2) Tokuyama Universal Bond II does not adhere to tooth structure immediately after applications of materials that contain:
 - hydrogen peroxide (oxydol),
 - sodium hypochlorite

We advise avoiding the use of these materials on the day of cementation, however these materials can be used for root canal treatment as the inhibition caused by these materials usually evaporates within five days.

■ **STORAGE**

- 1) Store Tokuyama Universal Bond II at temperatures between 0 to 25°C (32 to 77°F).
- 2) Keep away from heat, direct sunlight, sparks and open flames.
- 3) DO NOT use after the indicated expiration date on the bottle/package.

■ **DISPOSAL**

Unused Tokuyama Universal Bond II is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

■ **CLINICAL APPLICATION FOR DIRECT RESTORATIONS AND INTRAORAL REPAIR OF RESTORATIONS WITH COMPOSITE RESIN**

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Isolation

The rubber dam is the preferred method of isolation.

3. Cavity preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.

- If the adherent surface involves uncut enamel, apply etching agent to the uncut enamel. See section “Selective-enamel-etch technique”.

4. Drying

Dry the cavity by using a blotting technique or an air syringe technique.

- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.

- The substances listed below, which inhibit curing of Tokuyama Universal Bond II, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, or the application of phosphoric acid for 2 to 3 seconds before application:

- 1) Oil mist from handpiece,

2) Saliva, blood and exudates.

5. Etching of the tooth structure (optional)

(a.) Selective-enamel-etch technique

Apply phosphoric acid gel only to the enamel surrounding the margin of prepared cavity and leave phosphoric acid gel in place for 10-15 seconds. Rinse the etched surface thoroughly (at least 15 seconds) with water, and then dry with mild air.

(b.) Total etch technique

Apply phosphoric acid gel onto the prepared enamel first, and then the dentin. The etchant should be left to react on the enamel and dentin for 10-15 seconds. Rinse thoroughly (at least 15 seconds) with water and dry with mild air or with cotton pellets; do not overdry.

6. Preparing the adhesive surface during repair

6.1 Ceramics other than Porcelain, and Composite materials

Bevel the margin if necessary. Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

6.2 Metal restorations

Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

6.3 Porcelain

Bevel the margin if necessary. DO NOT roughen the surface; Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.

- When using air abrasion, use at least 30 to 50 micron grid of alumina powder.

- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking.

1) Oil mist from handpiece

2) Plaque, tartar, saliva, blood and exudate fluids

- If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

7. Dispensing of Tokuyama Universal Bond II

Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into the mixing well or disposable mixing well. Mix thoroughly with a disposable applicator until the mixed bonding agent turns green.

- Dispense Tokuyama Universal Bond II B IMMEDIATELY after dispensing A. Tokuyama

Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.

- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using disposable mixing well.
- Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.
- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the cap prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- DO NOT over tighten the bottle caps.

8. Application of Tokuyama Universal Bond II

Using the disposable applicator, apply the mixed Tokuyama Universal Bond II to the entire adherent surface(s). Take care not to omit any surfaces where the Tokuyama Universal Bond II should be applied.

- In the case of intraoral repair of restorations, if the adherent surface includes tooth structure and restorations (ceramics, porcelain, metals or composite materials), apply the mixed Tokuyama Universal Bond II to the entire adherent surface at once. No extra primers are needed.
- If the applied Tokuyama Universal Bond II pools on the cavity floor or a cavosurface angle and is too thick to air thin, remove the excess with a paper point or a new disposable applicator.
- Be careful not to allow excess Tokuyama Universal Bond II to flow into the subgingival. If Tokuyama Universal Bond II flows into the subgingival, thoroughly rinse and dry the area.
- If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse the cavity with water, dry and re-apply fresh Tokuyama Universal Bond II.

- DO NOT rinse the applied Tokuyama Universal Bond II with water except in unintentional contamination.
- DO NOT wait too long after the application of Tokuyama Universal Bond II. Air dry within 30 seconds after application, to ensure the appropriate film thickness.

9. Air Dry

Using an oil-free air/water syringe, apply weak air continuously to the Tokuyama Universal Bond II surface until the runny Tokuyama Universal Bond II stays in the same position without any movement. Finish by using mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.

- If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
- Do not rinse the applied Tokuyama Universal Bond II with water except in unintentional contamination.
- If accidental spattering occurs, it may cause the tissue to whiten or possible allergic reaction.

10. Composite resin

Restore with a light-cured, self-cured or dual-cured composite resin according to its manufacturer's instructions. The overfilled composite resin should be thoroughly finished and polished.

- The Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.

■ CLINICAL APPLICATION FOR CEMENTATION OF INDIRECT RESTORATIONS

1. Preparing the adhesive surface of the restoration

1.1 Ceramics other than Porcelain, and Composite materials

Roughen the interior of the restoration by sandblasting (0.1 to 0.2MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion.

Thoroughly rinse and dry the surface.

1.2 Metal restorations

Roughen the area by sandblasting (0.3 to 0.5MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

1.3 Porcelain

DO NOT roughen the surface; Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.

- When using the sandblaster or air abrasion, use at least 30 to 50 micron grid of alumina powder.
- If unable to roughen or polish the adhesive surface, etch the area for adhesion with phosphoric acid then thoroughly rinse and dry.
- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking. Clean the restoration after adjustment is complete.
 - 1) Gypsum and investment
 - 2) Silicone oil from fit check material
 - 3) Oil mist from handpiece
 - 4) Plaque, tartar, saliva, blood and exudate fluids
 - 5) Temporary restoration and temporary cement

If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

2. Dispensing of Tokuyama Universal Bond II

Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into the mixing well or disposable mixing well. Mix thoroughly with a disposable applicator until the mixed bonding agent turns green.

- Dispense Tokuyama Universal Bond II B IMMEDIATELY after dispensing A. Tokuyama Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.
- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using disposable mixing well.
- Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.
- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the

cap prior to closing.

- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- Do not over tighten the bottle caps.

3. Application of Tokuyama Universal Bond II

Using the disposable applicator, apply the mixed Tokuyama Universal Bond II on the surface to be bonded.

- Ensure Tokuyama Universal Bond II covers all surfaces where the cement will be applied.
- DO NOT over apply. One application of Tokuyama Universal Bond II is sufficient. Too thick of a layer may compromise the fit of the restoration.
- DO NOT allow saliva, blood, oil, water etc. to come in contact with the surface to be bonded. If the treated surface is contaminated, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
- DO NOT wait too long after the application of Tokuyama Universal Bond II. Air dry within 30 seconds after application, to ensure the appropriate film thickness.

4. Air dry

Using an oil-free air/water syringe, apply mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.

5. Cementing

Cement with a resin cement according to its manufacturer's instructions.

- DO NOT apply resin cement paste directly to the prepared tooth or the cavity treated with Tokuyama Universal Bond II to avoid improper seating of the restoration.
- Curing may be accelerated depending upon the resin cement used. To keep the appropriate film thickness, seat the restoration as soon as possible after the application of the resin cement (within 1 minute after dispensing), if not there is a possibility that the restoration will not seat properly with cementation.
- The Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.

CLINICAL APPLICATION FOR BONDING OF ACRYLIC RESIN TO METAL BASE, CLASP OR ATTACHMENT

1. Preparing the adhesive surface

Roughen the area by sandblasting (0.1 to 0.2MPa), or by grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

2. Dispensing, Application and Air dry

Use the same procedures as described in CLINICAL APPLICATION FOR CEMENTATION OF INDIRECT RESTORATIONS.

3. Application, curing, trimming and finishing of acrylic resin

Proceed with the procedure according to the manufacturer's instructions for the use of polymerizable acrylic resin.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.

Specifications are subject to change without notice. When the product specifications change, the instructions and precautions may also change.

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