**ESTESIL H<sub>2</sub>TOP®**

is a line of addition reaction impression materials in polyvinylsiloxane with highly wettable surfactant hydrophilic fillers. It is available in different formulations and viscosities to meet the requirements of the main impression taking techniques.

**WARNING: For dentistry usage only.**

**We suggest to make the plaster model at least 2 hours after taking the impression, in order to allow a perfect elastic recovery.**

**COMPOSITION**

Base: mixture of polyvinylsiloxane, polyvinylsiloxane terminal - H, inert fillers and color pigments.  
Catalyst: mixture of polyvinylsiloxane, inert fillers and patented platinum complex.

**TECHNICAL DATA**

CODE TYPE	ISO 4823 CLASSIFICATION							
	B PUTTY SOFT	Δ MONOPHASE	I MONOIMPLANT	E HEAVY	R REGULAR FLOWABILITY	F HIGH FLOWABILITY	H EXTRA HIGH FLOWABILITY	
PROPERTIES	Ratio	1:1	5:1	5:1	1:1	1:1	1:1	1:1
Colour	Yellow	Acid Green	Purple	Ochre	Grey	Blue	Pink	
Maximum Working Time	2 min 0 sec	2 min 30 sec	2 min 15 sec	2 min 30 sec	2 min 15 sec	2 min 15 sec	2 min 15 sec	
Total Setting Time	*3 min 30 sec	*4 min 0 sec	*4 min 0 sec	*4 min 0 sec	*4 min 0 sec	*4 min 0 sec	*4 min 0 sec	
Classification	ISO 4823	Type 0 (ADA)	Type 2 Medium	Type 2 Medium	Type 1 High	Type 2 Medium	Type 3 Light	Type 3 Light
Hardness SHORE A	Point after 24h ASTM D2240	65	60	58	62	55	53	50
Dimensional Stability	% ISO 4823	< 0,2	< 0,2	< 0,2	< 0,2	< 0,2	< 0,2	< 0,2
Elastic Recovery	% ISO 4823	99,8	99,8	99,7	99,7	99,8	99,7	99,8
Strain in Compression	% ISO 4823	1 - 3	3 - 5	3 - 5	3 - 5	3 - 5	3 - 5	3 - 5
Tear Strenght	N/m	NA	8	8	9	8	8	7
Compatibility with Paris plaster		Excellent	Excellent	Excellent	Excellent	Excellent	Excellent	Excellent

\* Recommended total setting time 4 min 30 Sec

**STORAGE**

An improper and inadequate storage may negatively affect the product characteristics causing its malfunctioning, or it may decrease its shelf life. The product packages must be stored away from heat sources, in a dry place with temperature ranging from 10°C/50°F to 24°C/75°F. Do not freeze the product packages. Do not use after the expiry date.

**SUGGESTION FOR THE IMPRESSION MATERIALS MATCHING**

TYPE	COLOR	CODE	SUGGESTIONS FOR USE
H Heavy	OCHRE	E	E+F: precision impression for supra-gingival or iuxta gingival preparations in the dual phase technique
RF Regular Flowability	GREY	R	B+R: precision impression for supragingival preparations when you need closer viscosities
HF High Flowability	BLUE	F	B+F or Δ+F or E+F: precision impression for supragingival or iuxta gingival preparations
EHF Extra High Flowability	PINK	H	B+H o Δ+H: precision impression for subgingival reading
PUTTY SOFT	YELLOW	B	Base for dual phase impression with two viscosities (putty-wash technique)
MONOPHASE DELTA	ACID GREEN	Δ	Base for simultaneous impression with two viscosities for: Dual Arch, Implant or Traditional. Recommended also for Inlays, Onlays
MONOIMPLANT	PURPLE	I	For monophasic impression (one step-single viscosity) for: Implant or Traditional

**ESTESIL H<sub>2</sub>TOP® Putty**

This is a addition reaction paste for manual mixing for impression taking in polyvinylsiloxane with highly wettable surfactant hydrophilic fillers, high viscosity and with mechanical properties ideal for impression taking in the dual phase two viscosity technique (putty-wash technique).

**INDICATIONS**

As a base in the dual phase two viscosity technique, of impressions for: crown, bridges, inlays, onlays and veneers.

**RECOMMENDATIONS TO USE ESTESIL H<sub>2</sub>TOP® Putty:**

Take out the same amount of base and catalyst using the special measures with the corresponding colors and mix with your fingers (approx 30 seconds) until the color looks homogeneous. Fill the tray immediately.

**WARNINGS**

- Talc powdered gloves or containing sulphur may interfere with the polymerization of A-silicones (or polyvinylsiloxanes). If you have used this type of gloves wash your hands with care.
- Do not exceed in the use of surfactants because they may inhibit polymerization.
- Astringents containing aluminum salts may inhibit the polymerization of A-silicones (or polyvinylsiloxanes).
- Do not exchange the components in the packaging and the measures to avoid affecting the product's properties.
- Do not use simultaneously and in mono or dual phase techniques other materials like: polyethers, polysulphurs and condensation silicones.

**ESTESIL H<sub>2</sub>TOP® Putty** can be used with most trays on the market with adequate retention perforations or retentions features. It is also possible to use a silicone adhesive (for example Tokuyama Adhesive), to enhance adhesion of the material to the stainless steel or plastic tray.

**IMPRESSION TAKING**

Dual phase two viscosities technique: fill the tray previously selected for the impression and insert it in the oral cavity in the time indicated in the table. At the end of the setting time indicated remove the tray from the oral cavity, checking the impression in the mouth (on the outer edge) to make sure that the curing and hardening of the material are concluded.

**IMPORTANT**

To achieve a better cohesion with the correction material do not wait more than 30' from the first impression taking.

**SUGGESTIONS FOR THE IMPRESSION REMOVAL**

Remove the tray from the oral cavity following the tooth axis; avoid to remove it on a slanted or tilted line because this may strain the material or cause imperfections.

**SUGGESTION FOR DISINFECTION**

To minimize the bacterial charge spray or immerse the impression in special solutions for addition reaction silicones, wait and observe the time recommended by the disinfectant's manufacturer. To minimize environmental impact, Tokuyama Dental recommends using a water-based solution disinfectants.

Remove the impression from the disinfectant solution, rinse thoroughly to remove all disinfectant residues, let dry for an adequate time (at least the same time of the disinfection time) and store properly at room temperature (approx. 25°C).

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**ESTESIL H<sub>2</sub>TOP® MONOPHASE Δ AND MONOIMPLANT**

**ESTESIL H<sub>2</sub>TOP® 5:1** is an addition reaction impression material in polyvinylsiloxane with highly wetttable surfactant hydrophilic fillers in the version with cartridge for mechanic mixing. Usable with the main impression material mixing machines, in a 5:1 ratio.

The indications are the following:

**INDICATIONS**

As a base in the simultaneous two viscosity technique of impressions for: crown, bridges, inlays, onlays and veneers. Indicated for the dual arch technique.

In the monophasic single viscosity of impressions for: positioning of pickups of implants, antagonists, position and functional impressions.

**PREPARATION**

Using the Mix Machine 5:1 cartridge proceed as follows:

1. Activate the cartridge removing the protection cap.
2. Insert the mixing device on the cartridge.
3. Fix the mixing device with the locking ring to the cartridge.
4. Insert the assembly in the special space in the mixing machine (follow the manufacturer's instructions on the use and position of the cartridges in the mixing machine).
5. Dispense about 3 cm. of material, and eliminate it, before placing it on the tray.
6. Do not remove the mixing device from the cartridge once it has been used (it acts as a protection plug).

**TRAY**

For this type of technique it is possible to use any type of non perforated trays with adequate retention features, individual or semi-individual trays and for the dual arch technique.

**ESTESIL H<sub>2</sub>TOP® HEAVY, REGULAR, HIGH AND EXTRA HIGH FLOWABILITY**

Are addition reaction impression materials in polyvinylsiloxane with highly wetttable surfactant hydrophilic fillers used as precision correctors. Available in different formulations to meet the requirements of the main impression taking technique. They are supplied in the cartridge version to be used with mixing guns in a 1:1 ratio.

Using the Gun Mix 1:1 cartridge proceed in the following way:

1. Insert the cartridge into the mixing gun.
2. Activate the cartridge removing the protection cap.
3. Insert the mixing cannula on the cartridge.
4. Dispense about 2/3 cm. of material, and eliminate it, before placing it in place directly with the special tips, or injecting it with special syringes used to achieve a precise application into the area where a more detailed impression taking is necessary.
5. Do not remove the mixing cannula once it has been used (it acts as a protection plug).

**GENERAL SAFETY INFORMATION**

The formulations of the impression materials **ESTESIL H<sub>2</sub>TOP®** contain curing monomers and surfactants which may irritate the skin, the eyes and the oral mucosa. They may cause allergic contact dermatitis in particularly susceptible subjects.

**CONTRAINDICATIONS**

1. Avoid contact with the eyes to prevent potential irritations and /or accidental damage to the cornea. In case of accidental contact with the eyes rinse with plenty of water and call a doctor.
2. Avoid contact with the skin to prevent potential irritations or cause any allergic reactions (in susceptible subjects). In case of contact with the skin remove the material and rinse with plenty of water and neutral soap. In case of sensitization or reddish rash, stop using the material and seek medical assistance.
3. Do not ingest the product. In case of accidental ingestion drink plenty of water. In general the material is not dangerous if swallowed in small quantities. Seek medical advice in case of symptoms of a different nature.

**PRECAUTIONS**

1. **ESTESIL H<sub>2</sub>TOP®** is intended only for dentistry, following the specific "instructions for use". Any different and/or imprudent use or and/or use not compliant with the specifications expressly provided is to be considered the dentist's sole discretion and responsibility.
2. Do not use dental materials containing hydrogen peroxide in combination with **ESTESIL H<sub>2</sub>TOP®** because they may prevent its correct hardening.
3. Thanks to its high thixotropy and formulation, the material should be dispensed easily. DO NOT apply too much pressure to avoid affecting the extrusion of the material or breaking the packaging.

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**ESTESIL H<sub>2</sub>TOP® HEAVY, REGULAR, HIGH ET EXTRA HIGH FLOWABILITY**

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For Tokuyama Dental Corporation
Tokuyama Dental Italy S.r.l.
Pour Tokuyama Dental Corporation
Per Tokuyama Dental Corporation